

## mediNiK® – An innovative hydrogel for removal of kidney stone fragments following lithotripsy


### WHY ARE SMALL RESIDUAL FRAGMENTS (RF) AND A STONE-FREE RATE OF 100% RELEVANT?

- **Nearly half** of all stone fragments after lithotripsy are < 1 mm and therefore **cannot be grasped**<sup>1</sup>
- **In 60%** of patients with RFs < 1 mm there was **no complete passage** of the RFs after 2 years and **in 18%** an **increase in size of the RFs**<sup>2</sup>
- **20%** of patients with RFs < 4 mm require an **intervention** within 20 months and there is **no significant difference in the disease progression rate** between patients with small (< 4mm) and large (> 4 mm) RFs<sup>3</sup>

Small fragments should be removed completely and a stone-free rate of 100% should be reached to reduce the rate of recurrence. To date these have been difficult to grasp.

### How does mediNiK® support the goal of stone-free status?

#### STONE REMOVAL RATES WITH mediNiK®




These results are from a non-published clinical investigation study. The study is an open-label, prospective, randomized multicenter study. It was conducted in the period from 09-2021 to 09-2022 in 5 investigation sites in Germany (Asklepios Klinik Barmbek (Hamburg), Universitätsklinikum Münster (Muenster), Krankenhaus der Barmherzigen Brüder (Trier), Klinikum rechts der Isar (Munich), and Marienhaus Klinikum (Ahrweiler)). 59 patients with kidney stones > 5mm were either treated with URS with laser lithotripsy and stone fragmentation (dusting and pop-dusting) alone or additionally with mediNiK®.<sup>4</sup>

With mediNiK® fragments < 1 mm can be grasped and removed more effectively compared with URS alone.<sup>4</sup>

# Advantages of mediNiK® compared to blood-clot method

## COMPARISON OF RETRIEVAL TIMES OF FRAGMENTS ≤ 1 MM\*\*



\*\* An ex vivo pig kidney model was used in which 30 kidney stone fragments ≤ 1 mm were introduced. Not all fragments could be retrieved in 40% of the experiments in the URS alone group and they were stopped after 45 min.<sup>5</sup>

With mediNiK® fragments < 1 mm can be retrieved significantly faster.<sup>5</sup>

## COMPARISON OF BLOOD-CLOT METHOD WITH mediNiK®<sup>6</sup>

	Blood-Clot methode	mediNiK®
Effectiveness of retrieval	✓	✓
Rapid clot/gel formation	—	✓
Good vision	—	✓
Removal of residual clot/gel by diuresis	—	✓

mediNiK® is superior to the blood-clot method.



Small fragments (< 1 mm) can be grasped.<sup>4</sup>



Small fragments are retrieved efficiently.<sup>4</sup>



Superior to the blood-clot method.<sup>5,6</sup>



### SOURCES —

- Reddy, N.K., et al., Size Distribution of Fragments by High-power Holmium Laser Lithotripsy in MiniPCNL with Suction. *Curr Urol Rep.* 2021 Dec 16;22(12):64. doi: 10.1007/s11934-021-01072-8. PMID: 34913152.
- Kang, M., et al., Clearance rates of residual stone fragments and dusts after endoscopic lithotripsy procedures using a holmium laser: 2-year follow-up results. *World J Urol.* 2016. 34(11): p. 1591-1597.
- Brain E, Geraghty RM, Lovegrove CE, Yang B, Somani BK. Natural History of Post-Treatment Kidney Stone Fragments: A Systematic Review and Meta-Analysis. *J Urol.* 2021 Sep;206(3):526-538. doi: 10.1097/JU.0000000000001836. Epub 2021 Apr 27. PMID: 33904756.
- Netsch C. & Simovic S. Open label, Randomized, Multicentric Study to Evaluate Safety, Tolerability and

- Performance of mediNiK in Comparison to Standard of Care in Removal of Kidney Stones (2023); [not published clinical investigation report]; investigation sites: Asklepios Klinik Barmbek (Hamburg, Germany), Universitätsklinikum Münster (Muenster, Germany), Krankenhaus der Barmherzigen Brüder (Trier, Germany), Klinikum rechts der Isar (Munich, Germany), and Marienhaus Klinikum (Ahrweiler, Germany).
- Schoeb, D.S., et al., New for Old-Coagulum Lithotomy vs a Novel Bioadhesive for Complete Removal of Stone Fragments in a Comparative Study in an Ex Vivo Porcine Model. *J Endourol.* 2017. 31(6): p. 611-616.
- Straub, M, Stone Debris Management: Magnetic particles for extracting fragments, [conference contribution], Technology & Training in Endourology, Torino, Italy, 23-25. November 2022

